

Richie Brace Rx Form

HERSCO
ORTHODONTIC LABS

800.301.8275

39-28 Crescent Street
Long Island City, NY 11101

Doctor Name: _____

Account Number: _____

Address: _____

City: _____ State: _____ Zip: _____

Patient Name: _____ [] Right

Age: _____ Ht: _____ Wt: _____ Shoe Type/Size: _____

Shoes Enclosed: Yes: No: Male: Female: [] Left

DIAGNOSIS

Describe patient condition

RICHIE BRACE® Prescription Please mark medial, lateral malleoli and accommodations on cast

Richie Brace® Prescription

- Richie Brace® Standard: -full flexion ankle pivot
- Richie Brace® Restricted Ankle Pivot: -limits ankle motion, yet allows smooth contact phase of gait
• Indications: DJD ankle & STJ, dropfoot, tarsal coalition, mild Charcot, lateral ankle instability, peroneal tendinopathy.

Color Options

- Black (standard)
- White
- Beige

- Richie Brace® Dynamic Assist: -full flexion pivot with spring hinges for dorsiflexion assist
• Patient requirements: 1: Dropfoot 2: Ankle dorsiflexion to at least 90° to leg 3: Stable knee ~ (must have all 3)
- Little Richie Brace® - Pediatric application for shoe size 4 and under
- Richie Soccer Brace® - Includes shin guard
- Richie Brace Ultra™ Modification (provide height/weight above)

- Richie Brace® Solid AFO: -Traditional full leg posterior shell w/ balanced functional orthotic footplate
• Indications: Dropfoot with unstable knee, dropfoot with spasticity, Charcot Arthropathy
STS BERMUDA CASTING SOCK REQUIRED

Color Options

- Tan
- Chocolate

- Richie California® [] 7" [] 9" (standard)
- Richie Gauntlet® [] 7" (standard) [] 9"
- BOTH GAUNTLET AND CALIFORNIA REQUIRE THE STS MID LEG SOCK / Has a medial arch suspender unless specified other wise

RICHIE BRACE® MODIFICATIONS - Note: Non-standard brace modifications may have extra charges - see pricing sheet

- Medial Arch Suspender -Adjustable lifting strap under talo-navicular joint for severe PTTD
- Lateral Arch Suspender -Adjustable lifting strap under calcaneal-cuboid joint for peroneal tendinopathy and severe lateral ankle instability
- Posterior Upright Connector -Connects uprights to stiffen brace (Arch suspenders require either a restricted ankle pivot or a posterior upright connector)

Top Cover	Length	Heel Cup	Medial Heel Skive
<input type="checkbox"/> EVA (standard)	<input type="checkbox"/> to Mets (standard)	<input type="checkbox"/> 10mm	For severe pronation control
<input type="checkbox"/> Terryco	<input type="checkbox"/> to Sulcus	<input type="checkbox"/> 14mm	<input type="checkbox"/> 2mm
<input type="checkbox"/> Diabetic (Plastizote/Poron)	<input type="checkbox"/> to Toes	<input type="checkbox"/> 18mm	<input type="checkbox"/> 4mm
	<input type="checkbox"/> 1/8" poron cushion on extension	<input type="checkbox"/> 35mm (standard)	<input type="checkbox"/> 6mm
Footplate Accommodation (please mark on cast)	Footplate Modification	Foot Plate Thickness	Extrinsic Posting - Rearfoot
<input type="checkbox"/> Navicular	<input type="checkbox"/> Medial Arch Flange	<input type="checkbox"/> 3.0mm (standard < 200 lbs)	<input type="checkbox"/> Heel Stabilizer Bar (standard)
<input type="checkbox"/> Styloid 5th Met	<input type="checkbox"/> Fascia Band	<input type="checkbox"/> 4.0mm (standard > 200 lbs)	<input type="checkbox"/> Rearfoot Post _____° Varus
<input type="checkbox"/> Other	<input type="checkbox"/> Lateral Flange	<input type="checkbox"/> 5.0mm	_____° Valgus
			<input type="checkbox"/> Heel Lift (Requires rearfoot post)
			_____ (inches)
Other Modifications		Limb Uprights	
<input type="checkbox"/> Crepe Plantar Arch Fill		<input type="checkbox"/> Align perpendicular to foot plate (standard)	
<input type="checkbox"/> Sulcus Wedge _____° Varus _____° Valgus		<input type="checkbox"/> Align 10° inverted to foot plate (<10% tibial varum)	

Special Instructions:

Accommodation location(s):
(mark on illustration and on cast)

